



Chamber Well Insert

Protects chamber well, easily washable with decontaminate

[Get a Quote](#)

Category: [Dose Calibrator Accessories](#)

Description

Description

Every Capintec dose calibrator in the CRC family already comes with this protective chamber well insert. If it becomes contaminated, you can easily remove it for washing with a decontaminate. Or, you can buy an extra chamber well insert to use until the other one decays to background levels.

The chamber well insert is made of clear, tough Plexiglas for improved durability, reduced breakage, and longer life.

Item Numbers

7300-2004

Liner

[Get a Quote](#)

Related products




Ultra Tech T-Mac System



CAP-Lift



CAP-MAC-S Moly Assay Canister For Syringes



CRP-200 Dose Tickets And Labels

