

Rad Scooter

Home / SPECT Imaging Products / Rad Scooter



Rad Scooter

Radscooter™ was designed to provide an alternative for transferring your source to the camera with minimal exposure without moving a 60 plus pound case or carrying a 30 pound Radshield™ around your Nuclear Medicine Department or Hospital. For storage purposes it requires no more area than a typical hard case and is available for BM01, BM02, and BM04 model sources only.

[Contact Your Distributor](#)

Share This Product



Physical Dimensions (Maximum):

7 inches wide (18 cm) x 30.2 inches long (77cm) x 32.5 inches high (83cm)

Standard Reference Information

- Hazmat Transportation Training Modules
- Dangerous Goods Declaration Template (PDF)
- Dangerous Goods Declaration Example (PDF)
- D.O.T. Flowchart


Related products



Dose Calibrator Sources "E" vial Products



Rod Sources



Spot Markers



Cobalt 57 Line Sources

Find your local distributor today!

[FIND DISTRIBUTOR](#)



RadQual, LLC

4137 Commerce Circle, Idaho Falls, ID 83401

Regulatory & Manufacturing:

p: (208) 524-5300

Sales & Product Distribution:

p: (208) 524-5300

e: sales@radqual.com

Quick Links



[Shipping & Returns](#)

[Contact Us](#)

[FAQs](#)

[Distributor Login](#)