Laminate disk sources are designed for use as reference sources in situations where less attenuation of the source material is required. Each laminate disk source is constructed using a high-strength laminate material and is available in two diameters:

- 0.2 x 1.0 inches (0.51 x 25.4 mm)
- 0.2 x 1.75 inches (0.51 x 44.5 mm)

The maximum deviation of the delivered activity from the nominal values listed is ±20%. A secondary calibration (against a NIST Traceable source) is available for an additional cost. This will provide a calibrated activity with a maximum deviation of ±5% from the labeled activity. Beta sources are not available for calibration.

Activities will not exceed the U.S. NRC Exempt Quantity limit.