

SPOT MARKERS

ERS spot markers are used for a range of purposes including patient orientation during a camera study, gamma camera calibration and quality control, gamma probe and NaI detectors calibration and performance check.

Epoxy marked active area is located in the center of the spot marker with a spherical diameter of 1 mm (0.039") and 1.8 mm (0.071"). Overall dimensions of ERS spot markers are 25.4 mm x 6.35 mm (1" x 0.25")



PRODUCT CODE	PRODUCT TYPE	ISOTOPE	ACTIVITY	
			SI	NON SI
SPT3180	Spot Marker	Na-22	0.37 MBq	10 uCi
SPT3181	Spot Marker	Na-22	3.7 MBq	100 uCi
SPT3980	Spot Marker	Gd-153	3.7 MBq	100 uCi
SPT3980-Set of 6	Spot Marker	Gd-153	3.7 MBq	100 uCi
SPT8080-100u	Spot Marker	Ge-68	3.7 MBq	100 uCi
SPT8080-1mCi	Spot Marker	Ge-68	37 MBq	1 mCi
SPT8080-200u	Spot Marker	Ge-68	7.4 MBq	200 uCi
SPT8080-250u	Spot Marker	Ge-68	9.25 MBq	250 uCi
SPT8080-25u	Spot Marker	Ge-68	0.925 MBq	up to 25 uCi
SPT8080-50u	Spot Marker	Ge-68	1.85 MBq	50 uCi
SPT8180-25u	Spot Marker	Na-22	0.925 MBq	up to 25 uCi
SPT8280-25u	Spot Marker	Co-57	0.925 MBq	up to 25 uCi
SPT8380-25u	Spot Marker	Ba-133	0.925 MBq	up to 25 uCi
SPT8480-25u	Spot Marker	Cs-137	0.925 MBq	up to 25 uCi
SPT8580-25u	Spot Marker	Co-60	0.925 MBq	up to 25 uCi
SPT8880-25u	Spot Marker	Eu-152	0.925 MBq	up to 25 uCi
SPT8180-50u	Spot Marker	Na-22	1.85 MBq	50 uCi
SPT8180-100u	Spot Marker	Na-22	3.7 MBq	100 uCi
SPT8180-200u	Spot Marker	Na-22	7.4 MBq	200 uCi
SPT8180-250u	Spot Marker	Na-22	9.25 MBq	250 uCi
SPT8180-1mCi	Spot Marker	Na-22	37 MBq	1 mCi
SPT8280-50u	Spot Marker	Co-57	1.85 MBq	50 uCi
SPT8280-100u	Spot Marker	Co-57	3.7 MBq	100 uCi
SPT8280-200u	Spot Marker	Co-57	7.4 MBq	200 uCi
SPT8280-250u	Spot Marker	Co-57	9.25 MBq	250 uCi
SPT8280-1mCi	Spot Marker	Co-57	37 MBq	1 mCi
SPT8380-50u	Spot Marker	Ba-133	1.85 MBq	50 uCi
SPT8380-100u	Spot Marker	Ba-133	3.7 MBq	100 uCi
SPT8480-50u	Spot Marker	Cs-137	1.85 MBq	50 uCi
SPT8480-100u	Spot Marker	Cs-137	3.7 MBq	100 uCi
SPT8480-200u	Spot Marker	Cs-137	7.4 MBq	200 uCi
SPT8480-250u	Spot Marker	Cs-137	9.25 MBq	250 uCi
SPT8580-50u	Spot Marker	Co-60	1.85 MBq	50 uCi
SPT8880-50u	Spot Marker	Eu-152	1.85 MBq	50 uCi
SPT1290	Spot Marker	Co-57	117 MBq	3 mCi
SPT1291	Spot Marker	Co-57	1.85 MBq	50 uCi